

# PCR Protocols Methods and Applications

Edited by  
David H. Gelfand, John J. Sninsky, and  
Thomas J. White

The polymerase chain reaction (PCR) is a powerful new method with widespread applications in diagnosis. With over fifty chapters of this unique, comprehensive benchtop reference, complete range of PCR methods and Equipment, reagents, and supplies are

## Key features include:

- Comprehensive study of PCR protocols, from multiplex PCR methodology
- Equipment, reagents, and supplies
- Experimental set-up
- Contamination

© 1990 by Academic Press, a Harcourt Brace Jovanovich Company

PCR Protocols A Guide to Methods and Applications

# PCR Protocols

A Guide to Methods and Applications

Edited by

Michael A. Innis  
David H. Gelfand  
John J. Sninsky  
Thomas J. White

# PCR PROTOCOLS

---

A GUIDE TO METHODS AND APPLICATIONS

Edited by

Michael A. Innis,  
David H. Gelfand, John J. Sninsky  
*Cetus Corporation, Emeryville, California*

Thomas J. White  
*Hoffmann-La Roche, Inc., Emeryville, California*



ACADEMIC PRESS, INC.  
Harcourt Brace Jovanovich, Publishers  
San Diego New York Boston  
London Sydney Tokyo Toronto

08238080-050394

08238080-050394

This book is printed on acid-free paper. (∞)

Copyright © 1990 by Academic Press, Inc.  
All Rights Reserved.

No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage and retrieval system, without permission in writing from the publisher.

Academic Press, Inc., San Diego, California 92101  
United Kingdom Edition published by  
Academic Press Limited, 24-28 Oval Road, London NW1 7DX

Library of Congress Cataloging-in-Publication Data

PCR protocols.

Includes index.

1. Polymerase chain reaction. 2. Gene amplification.

I. Innis, Michael A. [DNLM: 1. DNA Polymerases.

2. Gene Amplification--methods. 3. Genetic Engineering  
--methods. 4. RNA Polymerases. QH 442 P3479]

QP606.D46P36 1989 574.87'328 89-6938

ISBN 0-12-372180-6 (alk. paper)

ISBN 0-12-372181-4 (pbk. : alk. paper)

Printed in the United States of America

91 92 9 8 7 6 5 4 3

CO

Contr:  
Prefac

Part  
BASI

1. C  
A

2. A  
B

3. A  
B

4. B  
C

5. B  
C

6. C  
D

7. D  
E

8. E  
F

## PROCEDURES TO MINIMIZE PCR-PRODUCT CARRY-OVER

Shirley Kwok

The ability of PCR to produce large numbers of copies of a sequence from minute quantities of DNA necessitates that extreme care be taken to avoid false positives. Although false positives can result from sample-to-sample contamination, a more serious source of false positives is the carry-over of DNA from a previous amplification of the same target. Because of the large numbers of copies of amplified sequences, carry-over of even minute quantities of a PCR sample can lead to serious contamination problems. The following is a list of procedures that will help to minimize the carry-over of amplified DNA.

### Physical Separation of Pre- and Post-PCR Amplifications

To prevent carry-over of amplified DNA sequences, reactions should be set up in a separate room or containment unit such as a biosafety cabinet. A separate set of supplies and pipetting devices should be

142

PCR Protocol: A Guide to Methods and Applications  
Copyright © 1990 by Academic Press, Inc. All rights of reproduction in any form reserved.

dedicated for the specific use of sets taken to insure that amplified DNA. Reagents and supplies should be taken and must never be taken from a being performed. Similarly, devices be taken into the containment area.

### Aliquot Reagents

All reagents used in the PCR must be stored in an area that is free of PCR oligonucleotides used for amplification purified in a PCR-product-free environment to minimize the number of able to record the lots of reagents used it can be more easily traced.

### Positive Displacement Pipettes

Contamination of pipetting devices often of samples. For example, the pipettes often contaminated with radioisotopes. To eliminate cross-contamination, positive-displacement pipettes are used. Those manufactured by Rainin (Micro and plungers. The units are complete.

### Meticulous Laboratory Technique

Although carry-over of amplified sequences is a factor of the false positive, cross-contamination can also be a factor. Consequently, only in setting up the amplification.

08238080, 050394

sample handling, from sample collection to sample extraction. The following are additional precautions that should be taken:

1. Change gloves frequently.
2. Quick spin tubes before opening them.
3. Uncap and close tubes carefully to prevent aerosols.
4. Minimize sample handling.
5. Add nonsample components (mineral oil, dNTPs, primers, buffer, and enzyme) to the amplification reactions before the addition of sample DNA. Cap each tube after the addition of DNA before proceeding to the next sample.

## Judicious Selection of Controls

First, for use as a positive control, select a sample that amplifies weakly but consistently. The use of strong positives will result in the unnecessary generation of a large amount of amplified sequences. If plasmid DNA containing the target sequence is used as a positive control, it should be substantially diluted. Depending on the detection system used, as few as 100 copies of target will suffice as a positive control. Second, use well-characterized negative controls. The extreme sensitivity of PCR may enable the detection of nucleic acid sequence from a sample that is negative by all other criteria. Third, include multiple reagent controls with each amplification. Because the presence of a small number of molecules of PCR product in the reagents may lead to sporadic positive results, it is important to perform multiple reagent controls. The reagent controls should contain all the necessary components for PCR but without the addition of template DNA. This system has proved to be extremely sensitive in detecting the presence of contaminants, as the absence of exogenous DNA enables the efficient amplification of just a few molecules of contaminating sequence.

Although amplified products are most problematic, other potential sources of contamination/carry-over need to be considered, especially when additional manipulations of the amplified DNA are performed. The cloning of amplified product is a case in point. Often, the amount of target generated from an amplification is insufficient for direct cloning and requires re-amplification of the target. To minimize re-amplification of nonspecific products, the band of interest

is first separated on a gel, excised, and a subsequent amplification. potentially result in cross-contamination of the product. Precipitation. For example, gel in 1 N HCl to depurinate any residual device should be used to cause the surfaces of UV transilluminated, a sheet of plastic wrap should the gel from the surface of the controls that have been amplified with the sample of interest should not be preparative gel.

The list below highlights other

1. Plasmid or phage DNA contain
2. Purified restriction fragment of
3. Dot blot apparatus
4. Microtome blades
5. Centrifuges
6. Speed Vacs/vacuum bottles
7. Dry ice/ethanol baths

Other sources of contamination: dures (the preparation of samples; similar care) will most certainly here will serve as a guide in implementing minimize if not eradicate carry-over

08238080 050394